



Better lives through livestock

Sequence Data Quality Control

Kennedy Mwangi

International Livestock Research Institute (ILRI)

Viral Pathogen Genome Sequencing and Bioinformatics Analysis Training Workshop
6th – 17th May, 2024

ILRI
INTERNATIONAL
LIVESTOCK RESEARCH
INSTITUTE



Overview

- Should be first step!
 - What your data look like
 - Uses tools such as FastQC & MultiQC
- Removes:
 - Low quality bases
 - Low complexity sequences
 - Adaptor sequences

What Reads Do You Get



Single End Run
(one fastq file)



Paired End Run
(two fastq files)



FastQ Format Data

```
@HWUSI-EAS611:34:6669YAAXX:1:1:5069:1159 1:N:0:
TCGATAATACCGTTTTTTTCCGTTTGATGTTGATACCAT
+
DF=DBD<BBFGGGGGGGGBD@GGGD4@CA3CGG>DDD:D,B
@HWUSI-EAS611:34:6669YAAXX:1:1:5243:1158 1:N:0:
TATCTGTAGATTTACAGACTCAAATGTAAATATGCAGAG
+
IIHIIHIIIIIIIIIIIIIIIIIIIIIIIIIIHIIIIHIIIII
@HWUSI-EAS611:34:6669YAAXX:1:1:5266:1162 1:N:0:
GGAGGAAGTATCACTTCCTTGCCTGCCTCCTCTGGGGCCT
+
:GBGGGGGGGGGGGGDGGDEDGGDGGGGGDHHDHGHGHGBGG:GG
```


Read Quality – Phred Score

A quality value Q is an integer representation of the probability p that the corresponding base call is incorrect.

$$Q = -10 \log_{10} P \quad \longrightarrow \quad P = 10^{-\frac{Q}{10}}$$

Phred Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10,000	99.99%
50	1 in 100,000	99.999%
60	1 in 1,000,000	99.9999%

https://en.wikipedia.org/wiki/Phred_quality_score



ASCII Encoding

- Each number is converted to one symbol:

40 : @	90 : Z	141 : a
41 : A	91 : [142 : b
42 : B	92 : \	143 : c
43 : C	93 :]	144 : d
44 : D	94 : ^	145 : e
45 : E	95 : _	146 : f
... : : : ...

ASCII Encoding: cont...

ASCII_BASE=33 Illumina, Ion Torrent, PacBio and Sanger

Q	P_error	ASCII	Q	P_error	ASCII	Q	P_error	ASCII	Q	P_error	ASCII
0	1.00000	33 !	11	0.07943	44 ,	22	0.00631	55 7	33	0.00050	66 B
1	0.79433	34 "	12	0.06310	45 -	23	0.00501	56 8	34	0.00040	67 C
2	0.63096	35 #	13	0.05012	46 .	24	0.00398	57 9	35	0.00032	68 D
3	0.50119	36 \$	14	0.03981	47 /	25	0.00316	58 :	36	0.00025	69 E
4	0.39811	37 %	15	0.03162	48 0	26	0.00251	59 ;	37	0.00020	70 F
5	0.31623	38 &	16	0.02512	49 1	27	0.00200	60 <	38	0.00016	71 G
6	0.25119	39 '	17	0.01995	50 2	28	0.00158	61 =	39	0.00013	72 H
7	0.19953	40 (18	0.01585	51 3	29	0.00126	62 >	40	0.00010	73 I
8	0.15849	41)	19	0.01259	52 4	30	0.00100	63 ?	41	0.00008	74 J
9	0.12589	42 *	20	0.01000	53 5	31	0.00079	64 @	42	0.00006	75 K
10	0.10000	43 +	21	0.00794	54 6	32	0.00063	65 A			

$$Q = -10 \log_{10} P \quad \longrightarrow \quad P = 10^{\frac{-Q}{10}}$$


Read Quality: FastQC

FastQC Report

Summary

- ✓ [Basic Statistics](#)
- ✓ [Per base sequence quality](#)
- ! [Per tile sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ✗ [Per base sequence content](#)
- ! [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ! [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ✓ [Overrepresented sequences](#)
- ✓ [Adapter Content](#)

✓ Per base sequence quality



Reads raw fastq files


Performs multiple checks

- Pass/warn/fail
- Compares to genomic library


HTML Report

Read Quality: MultiQC


- Aggregates QC information from multiple samples
- Large number of programs supported
- Combined HTML report



Base Call Qualities – Per Cycle




Base Call Qualities – Per Cycle



Clean-up options

Trimming 3' end:


- Remove adapter read through
- Remove poor quality bases



Some quality issues may need to also remove specific reads



Despite issues may still be good enough for what is needed e.g. mapping

Per-Read Quality




- Are all reads equally affected?
- Is there a subset of reads which are always poor whilst others are good?

Measuring Read-through Adapters




Adapter removal

Before Trimming



After Trimming



Library Dependent QC Metrics


Some QC metrics will be influenced by what you are sequencing



Concern or Expected?



- GC Content
- Base Composition
- Duplication

Library GC Content



- Generic summary of library composition at a read level
- Expect a normally distributed set of values centred on the overall GC content


GC Content: Cont...



Specific Contamination with single sequence or closely related sequences

Artificial sequences, ribosomal RNA, contaminants


Library Base Composition




- For every chemistry cycle we can look at the number of ATGC we call
- For Libraries with random start positions the composition should be the same for all cycles

Bias Composition Throughout

Wrong Sequence




Amplicon




Proportional biases of bases at specific positions: Very low diversity

Bias Composition Throughout Cont...



Consistent disproportional expression of bases

Duplication



- How frequently the exact same sequence appears in your library
- For WGS expect most sequences to be unique

Duplication: Cont...

If the exact same sequence appears more than once it could be...

Technical:

ATCCGAGCTATTCGGCGAGCTCGCC

ATCCGAGCTATTCGGCGAGCTCGCC

ATCCGAGCTATTCGGCGAGCTCGCC

- PCR duplicates

Coincidental:

ATCCGAGCTATTCGGCGAGCTCGCC

ATCCGAGCTATTCGGCGAGCTCGCC


ATCCGAGCTATTCGGCGAGCTCGCC

- Deep sequencing
- Highly present sequences
- Restricted diversity libraries

Duplication: Cont...


RNA-Seq

Percent of seqs remaining if deduplicated 42.16%



Amplicon

Percent of seqs remaining if deduplicated 4.12%



Overrepresented Sequences

- Extreme duplication
- The exact same sequence is a significant proportion of the whole library (which might not be duplicated overall)
 - Poly Sequences
 - Specific Sequences

Sequence content across all bases

%T
%C
%A
%G

Position in read (bp)



ILRI
INTERNATIONAL
LIVESTOCK RESEARCH
INSTITUTE

CGIAR

Overrepresented Specific Sequences

- Normally artificial sequences (primers, adapters, vectors etc)
- Can search a database of known sequences to find matches

Sequence	Count	Percentage	Possible Source
GATCGGAAGAGCACACGTCTGAACTCCAGTCACCTTGTAATCTCGTATGC	17957	0.14359551756800035	TruSeq Adapter, Index 12 (100% over 50bp)

Acknowledgments

